Direction for use

1. Information

Megafill MH ceram. Universal light-curing Anterior-Posterior
-Microhybridcomposite with high ceramic glass-filler content

Shades:

Incisal,

Enamel: A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; C3; C4; D2; D3; D4

Dentin: A1; A2; A3; A3.5; B2; B3; D3

2. Composition

Filler: The filler based of Siliciumdioxid and high content of a special microglass, which is X-ray opaque (Strontium - Aluminium - Bor - Silikatglass)

Organic matrix:

2,2-Bis-4(2hydroxy-3-methacryloxy-propyloxy)-phenyl-propan;

7,9,9-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan,

-1, 16-di-yl-bis-methacrylat;

3,6-Dioxaoctamethylendimethacrylat

3. Indications

For restorations of cavity classes I, II, III, IV, V, aesthetic restorations, diastema, enamel hypoplasia and discoloration.

4. Contraindication

Placement of Megafill MH ceram is contraindicated if the patient is known to be allergic of the ingredients in Megafill MH ceram.

5. Side Effects

Irritations resulting from direct contact with dentin in proximity to pulp cannot be ruled out. In such cases the dentin should be shielded with a calcium hydroxide compound and a subfilling (e. g. glass ionomer cement) Application of eugenol containing material, specifically base liners, is not indicated, it will impair the polymerization.

6. Processing and Placement

-Tooth cleaning

The tooth and the adjacent distal surfaces has to be cleaned with a bristle brush and prophylactic paste.

-Shade Selection

Shade matching should preferably be carried out in daylight using the MEGADENTA Shade Guide. After curing the shade will be identical with the shade guide. The shades are corresponding to the Vita-Lumin-Vacuum Shade Guide.

-Cavity Preparation

The cavity is prepared and dried (rubber dam) in accordance with the general rules of composite filling treatment, the smear layer has not to be removed and unprepared dentine is cleaned with a cleaning paste. The prepared cavity is now rinsed with water and dried with oil free air. Deeper lying areas of the cavity can be filled with a suitable underfilling material.

-Base Liner

To protect exposed dentine the application of a glass ionomer liner, phosphate cement or calcium hydroxide base liner (deep cavities) is recommended, thus is reducing the danger of secondary caries at the margins.

-Cavity Finishing

The enamel margins are carefully finished with superfine diamond or carbide burs.

-Acid Etching

The cleaned enamel edges are etched with C-Cid. C-Cid is applied with the supplied brush. After 30-45 s C-Cid is applied on the dentine allowed to take effect for another 10-15 s and then thoroughly rinsed with water. Excessive drying should be avoided in this case, as this would cause the exposed collagen fibre to collapse

-Application of the Primer

The Primer has now to be used after the recommendations of the Manufacture. Do not rinse off the Dentin primer!

The cavity is dried thoroughly with oil and water free air or until the dentine surface appears dull and dry. It is essential that the primed dentine and the etched enamel are dry and contaminant free for the bond application.

C-Bond is applied onto the primed Dentine in the dried cavity and the etched enamel, blown out thinly with oil free air and then cured for 10 s with a blue light source (e.g. Megalux CS).

Following the polymerization of C-Bond, a light-curing composite (e.g. Megafill MH ceram) can be applied.

-Application

The quantity required for the cavity is extruded from the syringe and placed with an instruments into the cavity. The syringe should be capped and turned back immediately after extrusion. After completed contouring the restoration is cured for 40s with Megalux CS unit or an other halogen light unit. The beam has to be directed to all surfaces of the restoration.

Restorations deeper than 2,5mm should be filled in layers, each layer has to be cured separately (multi-phase technique). The influence of oxygen in the air causes a thin non-polymerized layer at the surface called

"dispersion or smear layer". This layer should not to be removed as it will facilitate a solid chemical bond with the next layer of composite.

-Shaping and Finishing

The final forming and the removal of excess composite are carried out with rotating finishing instruments. For contouring and finishing fine and ultra fine finishing diamonds are used. For polishing flexible discs and polishing strips are suitable. All shaping and finishing have to be carried out under tooth preserving and with water-cooling. A sealing of the filling with C-Bond is recommendable.

7. Shelf life

3 years

8. Precautions and Storage

The etching contains 35% phosphoric acid. Contact with oral soft tissues, skin and eyes must be avoided.

Not to be stored above 25°C. Storage in refrigerator prolongs shelf life

9. Deliver units

Megafill MH ceram	single syringe	4,5 g
Megafill MH ceram	assortments with	3 syringes + accessories
		4 syringes + accessories
		6 syringes + accessories
		8 syringes + accessories

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Universal-Microhybrid-composite
with high ceramic glass filler
for Anterior-Posterior restaurations

MEGADENTA Dentalprodukte GmbH
D-01454 Radeberg, Germany
Telefon +49(0)3528 453-0 Fax +49(0)3528 453-21
Http://www.megadenta.de info@megadenta.de

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